

About Guardant360® Liquid CDx

Guardant360 Liquid CDx is a fast, FDA-approved comprehensive liquid biopsy for advanced solid tumors that provides results in 7 days*.

Test specifications

Sample input

Two 10mL Streck tubes of whole blood†

Storage and shipping

Do not freeze or refrigerate. Ship same or next day at room temperature

Turnaround time

7 calendar days from sample receipt to results*



*Median turnaround time from sample receipt to first results

†Minimum 5 mL in one tube required to process

Gene panel

All four major classes of alterations detected

SNVs and indels (741 genes)

ABCB1	ABRX	CARD11	CDKN2B	DDR1*	ELAVL2*	FAAP100	FH	GPATCH8*	IKBKE*	KRAS	MGA	NOTCH3	PDCD1	PPM1D
ABL1	AURKA	CASP8*	CDKN2C	DDR2	ELF3*	FAAP20	FLCN	GPC3*	IKZF1*	LATS1	MITF*	NOTCH4*	PDCD1LG2	PPP2CA*
ABL2*	AURKB	CASR	CEBPA	DDX17*	ELOC*	FAAP24*	FLT1	GREM1	IL1R1	LGR4	MKNK1*	NOVA1*	PDE7A	PPP2R1A
ABRAXAS1	AURKC	CAV1	CELF4*	DDX18*	EML4	FANCA‡	FLT3	GRIN2A*	IL2RA	LGR5	MLH1‡	NPM1*	PDGFRA	PPP2R2A
ACVR1	AXIN1	CBFB*	CEP295*	DDX27*	EMSY	FANCB*	FLT4*	GSK3B	IL2RB	LGR6	MLH3*	NPRL2	PDGFRB	PPP3CA
ACVR1B*	AXIN2	CBL	CFAP20*	DDX3X*	EP300	FANCC*	FOXA1*	GSTM1	IL2RG	LIG1*	MLST8	NPRL3	PDK1*	PPP6C*
ACVR2A*	AXL	CBLB*	CHD4	DDX41*	EPAS1*	FANCD2*	FOXO1*	GSTP1	IL7R	LIG4	MPL	NRAS	PDPK1	PRDM1
ADARB2	B2M*	CCAR1*	CHEK1	DEPDC5	EPCAM	FANCE*	FOXO1*	H3-4	INHBA	LMO1*	MRAS*	NRG1	PHF6*	PRES1
ADGRA2*	BABAM1	CCN6*	CHEK2‡	DEPTOR	EPHA3	FANCF	FOXO1*	H3F3A*	INPP4B	LRP1B	MRE11	NSD1*	PHLPP1	PRES2
ADGRG4*	BABAM2	CCNA2	CIC	DHX15	EPHA5	FANCG*	FRS2	HACD4*	INTS6L*	LRP2*	MSH2‡	NSD2	PHLPP2	PRKAR1A
AFDN*	BAP1	CCNB1	CLDN18	DHX16*	EPHA7*	FANCI	FUBP1*	HDAC2	IRF1	LRP5	MSH3*	NSD3	PHOX2B	PRKCI*
AGGF1*	BARD1‡	CCND1	CMTM4	DHX36*	EPHB1*	FANCL	FUBP3*	HDAC6	IRF2	LRP6	MSH6*	NSRP1*	PIAS4	PRKDC*
AIP	BCL2	CCND2	CMTM6	DHX9*	ERBB2	FANCM*	FUS*	HELQ	IRF4*	LTK	MTAP	NTHL1	PIK3C2B*	PRKN
AKT1	BCL2L1	CCND3	CNOT3*	DICER1	ERBB3	FAS*	FYN	HES1	IRS2*	LYN	MTHFR	NTRK1	PIK3CA	PRMT5
AKT1S1	BCL2L2*	CCNE1	CREBBP*	DIS3L2	ERBB4	FAT1*	FZD1	HEY1	JAK1	LZTR1	MTOR	NTRK2*	PIK3CB	PRPF40B*
AKT2	BCL6*	CCNE2	CRKL*	DLL4	ERCC1*	FBXW7	FZD10	HEYL	JAK2	MAD2L2	MUTYH‡	NTRK3*	PIK3CD*	PRPF4B*
AKT3	BCOR	CD274	CRTC1*	DNAJB1	ERCC2	FCGR2A	FZD2	HGF	JAK3	MALT1	MYB*	NUMA1	PIK3CG	PSENEN
ALB	BCORL1*	CD276	CSF1R	DNMT1*	ERCC3*	FCGR3A	FZD3	HNF1A	JUN	MAP2K1	MYC	NUMB	PIK3R1	PSMB10
ALK	BCR*	CD74*	CSF3R	DNMT3A	ERCC4*	FEN1	FZD4	HNRNPDL*	KAT6A*	MAP2K2	MYCL	NUP93*	PIK3R2*	PSMB8
ALOX12B	BIRC5*	CD79A	CTC1	DNMT3B	ERCC5	FGF1	FZD5	HOXB 13	KAT6B	MAP2K4	MYCN	NUTM1*	PIK3R3	PSMB9
ALOX15B	BLM	CD79B	CTCF	DOT1L*	ERCC6	FGF10*	FZD6	HRAS	KDM4A	MAP3K1	MYD88	P2RY8*	PIM1	PTCH1
ALOX5	BMPR1A‡	CDC27	CTLA4	DPYD	ERCC6L2	FGF12	FZD7	HSD3B1*	KDM5A*	MAP3K13	MYOD1	PABPC1*	PIN1	PTDSS1
AMER1	BRAF	CDC5L*	CTNNA1	DUSP4	ERCC8	FGF14*	FZD8	HSP90AA1*	KDM5B*	MAP4K3*	NAB2	PAK1	PKM*	PTEN‡
APC‡	BRCA1‡	CDC7	CTNNB1	DYNLL1*	EREG	FGF19	FZD9	ICOSLG	KDM5C	MAPK1	NBN	PAK3*	PLCG2	PTPN11
APEX1	BRCA2‡	CDC73	CUL3	DYRK2	ERF	FGF2	GAS6	ID3	KDM6A	MAPK3	NCOR1*	PALB2‡	PLEKHS1	PTPN2
APLN	BRCC3	CDH1‡	CUL4A	E2F3*	ERG	FGF23*	GATA1	IDH1	KDR	MAPKAP1	NCR1	PARG	PLRG1*	PTPRD*
AR	BRD2*	CDH6	CUX1*	ECT2L*	ERRF1*	FGF3	GATA2	IDH2	KEAP1	MARK2	NCR3	PARP1*	PMS1*	PTPRS
ARAF	BRD3*	CDK11A*	CWC22*	EFTUD2*	ESR1	FGF4*	GATA3	IDO1	KIN*	MAX‡	NEGR1	PARP2*	PMS2‡	PTPRT
ARFRP1*	BRD4*	CDK12	CXCR4	EGFR‡	ETS1	FGF5	GATA4	IFNG	KIT	MCL1	NELFE*	PAX3	POLA1	QKI*
ARHGAP35	BRIP1‡	CDK4	CYLD*	EIF1AX*	ETV1	FGF6*	GATA6	IFNGR1	KLF4	MDC1*	NF1	PAX5	POLD1	RAB35
ARID1A	BSG	CDK6	CYP17A1	EIF4A1	ETV4	FGF7*	GEN1	IFNGR2	KLHL6*	MDM2	NF2‡	PAX7	POLE	RAC1*
ARID1B	BTG1*	CDK7	CYP19A1	EIF4A2	ETV5*	FGF8	GID4*	IFNW1	KLHL9	MDM4	NFE2L2	PAX8	POLH	RAD18*
ARID2	BTG2*	CDK8*	CYP2C19	EIF4A3	ETV6	FGF9	GLI1	IGF1	KMT2A	MED12	NFKBIA*	PAXIP1	POLQ*	RAD21*
ASXL1	BTK	CDKN1A	CYP3A4*	EIF4B	EWSR1*	FGFR1	GNA11	IGF1R	KMT2B	MEF2B*	NHEJ1*	PBRM1*	POT1	RAD50
ATM‡	BUB1B*	CDKN1B	DAXX	EIF4E	EXO1*	FGFR2	GNA13*	IGF2	KMT2C	MEN1‡	NKX2-1	PCBP1	POU2F2	RAD51
ATMIN	C9orf78*	CDKN1C	DCUN1D1	EIF4E2	EZH1	FGFR3	GNAQ	IGF2BP3*	KMT2D	MERTK	NOTCH1	PCBP2*	PPARG	RAD51B
ATR	CALR	CDKN2A‡	DDIT3*	ELAVL1*	EZH2	FGFR4	GNAS	IGF2R	KNSTRN*	MET	NOTCH2	PCDH15*	PPIG*	RAD51C*

*Regions of interest covered by the panel; some exons are not included.

‡Reporting is enabled for pathogenic germline alterations only. Somatic alterations are not reported.

*Reporting is enabled for both pathogenic germline and somatic alterations.

Gene panel (Continued)

All four major classes of alterations detected

SNVs and indels (741 genes)							CNAs (2 Genes)	CNL (1 Gene)	Fusions/Rearrangements (9 Genes)	
RAD51D [‡]	ROS1	SF3B1	SOX9	TEK	TPMT	WEE1	ERBB2 MET	BRCA1	ALK RET ROS1 NRG1 NTRK1	NTRK2 NTRK3 FGFR2 FGFR3
RAD52*	RPA1*	SF3B3*	SPEN	TEN1	TRAF2	WRN				
RAD54L	RPS27A*	SH2D1A	SPOP	TENT5C	TRAF3	WT1*				
RAET1E	RPS6KA3	SHLD1	SRC	TERT	TRAF7	WWP1				
RAF1	RPS6KB1	SHLD2*	SRSF2	TET1	TRIM24*	XBP1*	XPO1			
RARA	RPS6KB2	SLC34A2	SRY	TET2	TRIP13	XPA				
RASA1*	RPTOR	SLFN11	SS18	TFE3	TSC1 [‡]	XPC*				
RB1*	RRAGC	SLIT2*	STAG2	TFRC	TSC2 [‡]	XPO1				
RBBP6*	RSPO1	SMAD2	STAT1	TGFBF1	TSHR*	XRCC1*	XRCC2 XRCC3 XRCC4* XRCC5* XRCC6*			
RBM10	RSPO2	SMAD3*	STAT3*	TGFBF2	TSHZ2*	XRCC2				
RBMX*	RSPO4	SMAD4 [‡]	STAT4*	THRAP3*	TYMP	XRCC3				
RECQL	RUNX1	SMARCA2	STK11 [‡]	TIA1*	TYMS	XRCC4*				
RECQL4	RUNX1T1*	SMARCA4	STK19*	TIPARP	TYRO3*	XRCC5*	YAP1*			
RET [‡]	RXRA	SMARCAL1	STK40	TMEM127	U2AF1	XRCC6*				
REV3L	RYBP	SMARCB1	STN1	TMPRSS2	UBE2T*	YAP1*				
RGS1	SAMHD1	SMARCD1	SUFU	TNFAIP3	UGT1A1*	YES1				
RHEB	SDC4	SMARCE1*	SYK*	TNFRSF14*	UIMC1	ZC3H13*	ULBP1 ULBP3 USP28 USP7 USP9X*			
RHOA	SDHA [‡]	SMC1A	SYNCRIP*	TNFRSF1A	ULBP1	ZC3H18*				
RHOB	SDHAF2	SMC3	TACSTD2	TNK2*	ULBP3	ZC3H4*				
RICTOR	SDHB [‡]	SMO	TAF1L	TNPO1*	USP28	ZMYM3*				
RIF1	SDHC [‡]	SNCAIP	TAP1	TOP1	USP7	ZNF217	ZNF703*			
RILPL1*	SDHD [‡]	SOC1	TAP2	TOP2A	USP9X*	ZNF703*				
RIT1	SEM1*	SOC3	TAPBP	TOPAZ1*	VEGFA*	ZNRF3				
RNASEH2B	SERPINB3	SOS1	TBC1D7	TP53 [‡]	VEGFB	ZRSR2*				
RNF43	SERPINB4	SOX10*	TBX3*	TP53BP1	VHL [‡]		VIRMA*			
ROBO1	SESN2	SOX17	TCERG1*	TP63*	VIRMA*					
ROBO2	SETD2	SOX2	TCF7L2	TP73*	WBP11*					

*Regions of interest covered by the panel; some exons are not included.

[‡]Reporting is enabled for pathogenic germline alterations only. Somatic alterations are not reported.

[‡]Reporting is enabled for both pathogenic germline and somatic alterations.

Performance specifications

Alterations	Analytical Sensitivity*	Allelic Fraction/ Copy Number [‡]		Analytical Specificity*	Threshold for Positivity
		Low Input	High Input		
SNVs	≥ 94.9%	≥ 1%	≥ 0.2%	≥ 99.8%	≥ 0.001% MAF
Indels	≥ 95.0%	≥ 0.9%	≥ 0.2%	≥ 99.9%	≥ 0.01% MAF
CNAs	≥ 91.6%	≥ 2.3-2.4 copies	≥ 2.3-2.4 copies	≥ 99.9%	≥ 2.18-2.2 copies
CNL	≥ 95.0%	22.7% [§]	-	≥ 97.5%	Tumor fraction ≥ 20%
Fusions/Rearrangements	≥ 94.6%	≥ 0.6-1.6%	≥ 0.2-0.4% [§]	≥ 99.7%	≥ 2 unique molecules

*Analytical Sensitivity defined as analytical concordance with an orthogonal NGS device for clinically significant variants

[‡]Demonstrated Allelic Fraction/Copy Number at 95% detection rate, limit of detection (LoD), at low and high input

*Analytical Specificity defined as per-variant NPA against an orthogonal NGS comparator

[§]LoD established and/or verified using clinical samples for all variant categories

[‡]Reported as TF

Low input is defined as 400-500 NSC. High input is defined as 30ng.

Guardant360® Liquid CDx Intended Use

Guardant360® Liquid CDx is a qualitative next generation sequencing-based *in vitro* diagnostic device that uses targeted high throughput hybridization-based capture technology for detection of single nucleotide variants (SNVs) and insertions and deletions (indels) in 741 genes, copy number amplifications (CNAs) in two genes, copy number loss (CNL) in one gene, rearrangements in nine genes. Guardant360 Liquid CDx utilizes circulating cell-free DNA (cfDNA) from plasma of peripheral whole blood collected in Streck Cell-Free DNA Blood Collection Tubes (BCTs). The test is intended to be used as a companion diagnostic to identify patients who may benefit from treatment with the therapies listed in Table 1-1 in accordance with the approved therapeutic product labeling. Additionally, the test is intended to provide tumor mutation profiling to be used by qualified health care professionals in accordance with professional guidelines in oncology for cancer patients with solid malignant neoplasm. The test is for use with patients previously diagnosed with cancer and in conjunction with other laboratory and clinical findings.

Table 1-1 Companion Diagnostic Indications:

Indication	Biomarker	Therapy
Breast cancer	<i>ESR1</i> missense mutations between codons 310-547	ORSERDU™ (elacestrant)
	<i>ESR1</i> , <i>E380</i> , <i>V422del</i> , <i>S463</i> , <i>L469</i> , <i>L536</i> , <i>Y537</i> , and <i>D538</i> mutations	INLURIYO™ (imlunestrant)
Colorectal cancer	<i>BRAF</i> V600E	BRAFTOVI® (encorafenib) in combination with ERBITUX® (cetuximab)
Non-small cell lung cancer (NSCLC)	<i>EGFR</i> exon 19 deletions, L858R, and T790M*	TAGRISSO® (osimertinib)
	<i>EGFR</i> exon 20 insertions	RYBREVENT® (amivantamab-vmjw)
	<i>ERBB2/HER2</i> activating mutations (SNVs and exon 20 insertions)	ENHERTU® (fam-trastuzumab deruxtecan-nxki)
	<i>KRAS</i> G12C	LUMAKRAS™ (sotorasib)

A negative result from a plasma specimen does not assure that the patient's tumor is negative for genomic findings. Patients who are negative for the biomarkers listed in **Table 1-1** should be reflexed to tissue biopsy testing for **Table 1-1** biomarkers using an FDA-approved tumor tissue test, if feasible. *The efficacy of TAGRISSO (osimertinib) has not been established in the *EGFR* T790M plasma-positive, tissue-negative or unknown population and clinical data for T790M plasma-positive patients are limited; therefore, testing using plasma specimens is most appropriate for consideration in patients from whom a tumor biopsy cannot be obtained. Genomic findings other than those listed in **Table 1-1** are not prescriptive or conclusive for labeled use of any specific therapeutic product. Guardant360 Liquid CDx is a single-site assay performed at Guardant Health, Inc.

Contraindications

There are no known contraindications.

Warnings and Precautions

- Alterations reported may include somatic (not inherited) or germline (inherited) alterations. The assay filters germline variants from reporting except for pathogenic alterations in 38 genes. However, if a reported alteration is suspected to be germline, confirmatory testing should be performed using a validated germline test, as this assay is not intended to replace comprehensive germline testing or to provide definitive information about cancer predisposition.
- Genomic findings from cfDNA may originate from circulating tumor DNA (ctDNA) fragments, germline alterations, or non-tumor somatic alterations, such as clonal hematopoiesis of indeterminate potential (CHIP).
- Allow the Streck blood collection tube to fill completely until blood stops flowing into the tube. Underfilling of tubes with less than 5 mL of blood (bottom of the label indicates 5 mL fill when tube is held vertically) may lead to incorrect analytical results or poor product performance. This tube has been designed to fill with 10 mL of blood.

Limitations

- For *in vitro* diagnostic use.
- For prescription use only. This test must be ordered by a qualified medical professional in accordance with clinical laboratory regulations.
- The test is not intended to be used for standalone diagnostic purposes.
- The test is intended to be performed on specific serial number-controlled instruments by Guardant Health, Inc.
- A negative result for any given variant does not preclude the presence of this variant in tumor tissue.
- The efficacy of TAGRISSO (osimertinib) has not been established in the *EGFR* T790M plasma-positive, tissue-negative or unknown population and clinical data for T790M plasma-positive patients are limited; therefore, testing using plasma specimens is most appropriate for consideration in patients from whom a tumor biopsy cannot be obtained.
- ORSERDU efficacy has not been established in patients with *ESR1* missense mutations < 0.03% MAF.
- INLURIYO efficacy has not been established in patients with *ESR1* V422 deletion < 0.19% MAF and in patients with *ESR1* SNVs < 0.03% MAF.
- BRAFTOVI + ERBITUX efficacy has not been established in patients with *BRAF* V600E mutation < 0.06% MAF.
- TAGRISSO efficacy has not been established in patients with *EGFR* exon 19 deletions < 0.07% MAF, in patients with *EGFR* L858R < 0.12% MAF, and in patients with *EGFR* T790M < 0.05% MAF.
- RYBREVENT efficacy has not been established in patients with *EGFR* exon 20 insertions < 0.05% MAF.
- LUMAKRAS efficacy has not been established in patients with *KRAS* G12C biomarkers < 0.05% MAF.
- ENHERTU efficacy has not been established in patients with *ERBB2* exon 20 insertions < 0.03% MAF and in patients with *ERBB2* SNVs < 0.18% MAF.

Limitations (Continued)

- Decisions on patient care and treatment must be based on the independent medical judgment of the treating physician, taking into consideration all applicable information concerning the patient's condition, such as patient and family history, physical examinations, information from other diagnostic tests, and patient preferences, in accordance with the standard of care.
- ctDNA shedding rate may be lower in patients with primary central nervous system (CNS) tumors.



As a professional service, Guardant360 Liquid CDx assesses 746 genes and generates a professional services report that includes MSI-High status, bTMB, Viral (HBV/EBV) detection, pharmacogenomics, HLA genotyping, and additional features powered by InfinityAI.[^]

[^]These additional features are reported by Guardant360 Liquid CDx as a professional service and have not been cleared or approved by the US FDA.